

TELEMEDICINE CART

SOLUTIONS

2025

 Medical Cart
COMPANY



GENERAL FEATURES



EASY TO CLEAN SURFACES

Clean surfaces without dirt traps.
One wipe cleaning.



QUICK SETUP

Available assembled or
flat-packed, easy to assemble.



8 YEARS WARRANTY

Eight-year warranty for
the care-free treatment.



ERGONOMIC

Designed for easy
and comfortable use.



SMALL FOOTPRINT

Developed to fit narrow doors,
ease of movement.



LOCKABLE WHEELS

Simple locking mechanism
operated by foot.

CORE PRODUCTS

3-0302

HEIGHT ADJUSTABLE TELEMEDICINE CART


- Electric lift column
- B6 base - 80kg max load capacity,
125mm wheels
- Camera mount
- Display mount
- Shelf
- Powered base
- Inner cabling



3-0312

FIX HEIGHT TELEMEDICINE CART

- Fix height column
- B6 base
 - 80kg max load capacity
- Camera mount
- Display mount
- Powered base
- Inner cabling



3-0313

LIGHTWEIGHT TELEMEDICINE CART

- Fix height column
- B3 base
 - 35kg max load capacity
- Camera mount
- Display mount
- Keyboard tray
- CPU mount
- PDU
- DIN rail



FEATURES

HANDLE

Comfortable handle embracing the cart from the back. Lift actuator is placed on its center.



POWERED BASE

Battery pack can be easily mounted from below.



POWER DATA

Battery technology	LiFePO4
Usable WH / Battery capacity	450 WH / 512 WH
Battery life	2000 to 5000 charge cycles
Maximum load	250 W
Safety options	Chrgae indicator LED; Digital feedback on display

CABLE MANAGEMENT

The rear cover of the inner column can be removed to access the inside of the cart.



SAFE ACCESS

Key lock on the back. Small colored handles to remove the door. One cable output on the bottom.

DIGITAL FEEDBACK PANEL

Placed on the column back.
Service related conditions shown on the display.




DISPLAY AND CAMERA MOUNT

Two axis camera mount with hidden cabling.
Quick release display mount with rotate, tilt and turn functions.


Display mount Datasheet

Scan the QR



ACCESSORIES

VESA MOUNT




VES 1/2/3

Datasheet
Scan the QR



TABLET MOUNT



TM2

Datasheet
Scan the QR



CABLE HOOK



CH 1/2/3/4

Datasheet
Scan the QR



STORAGE




ST 1/2/3/5/6/7/8

Datasheet
Scan the QR



SHELVES



S 1/2/3/4/6/8/10/11

Datasheet
Scan the QR



DRAWERS




D 1/2/3/4/5/6/7/8/9/10/11

Datasheet
Scan the QR



SCANNER HOLDER




SH 1/2/3/4

Datasheet
Scan the QR



PSU HOLDER



PSH 1/2/4/5/7/8/9

Datasheet
Scan the QR



POWER STRIP



PDU 3/4/5/8

Datasheet
Scan the QR



COMPUTER MOUNT



CPU 1/2/3/4

Datasheet
Scan the QR



DIN RAIL




DR3, DR4, DR8, DR9

Datasheet
Scan the QR



POLE CLAMP



DRC 1/2

Datasheet
Scan the QR



ACCESSORIES

MOUSE HOLDER




MH1

Datasheet
Scan the QR



HUMIDIFIER MOUNT



AM1

Datasheet
Scan the QR



POAG STUD



POW-A2-15/25

Datasheet
Scan the QR



STANDARD COLORS

RAL 9005BM
Azzurro 3602
TIGER Drylac®
68/40413

RAL-9016
TRAFFIC WHITE

RAL 200-60-30
AQUA MARINE
TIGER Drylac®
14/40079

CUSTOM COLORS

We offer the choice of custom colors for the aluminium components from **Tiger 3D metallic range**.

Visit the Tiger webshop and select the color you desire:

<https://www.tiger-coatings.com>



Medical Cart C O M P A N Y



Medical Cart Company Kft.

2051 Hungary

Biatorbágy, Budai út 16.

Phone: +36 23 920 790

Email: info@medicalcartcompany.hu

www.medicalcartcompany.hu