



INTRAORAL CART


SOLUTIONS 2024



BEAVER 1

BV1


HEIGHT ADJUSTABLE
WORK STATION



BEAVER 1 PLUS

BV1-6-10


HEIGHT ADJUSTABLE
WORK STATION
WITH COMPUTER MOUNT



BEAVER 2

BV2


AIO COMPUTER
WORK STAND



BEAVER 2 PLUS

BV2-6-4-10


WORK STATION
WITH COMPUTER MOUNT



BEAVER 3

BV3-6-1


ALL-IN-ONE COMPUTER
WORK STATION



BEAVER 4

BV4-4-3


LAPTOP CART



BEAVER 5

2311-0618


SMALL DEVICE CART



BEAVER 6

BV6-4-9

PIEZOSURGERY CART



BEAVER 1

BV1

HEIGHT ADJUSTABLE WORKSTATION



Ergonomic



Easy to clean



Small footprint



8 years warranty




Quick setup



Locking wheels


Datasheet:
Click or Scan



FEATURES

HEIGHT ADJUSTMENT

Easy adjustment with gas spring.
Pull the lever upward to release
the lock and set it at the desired
height.



INTRAORAL SCANNER HOLDER (OPTIONAL)

We installed a holder
on the display mount for
your IntraOral scanner.

DISPLAY MOUNT

Tiltable and rotatable
VESA mount.
For displays up to 24".



BEAVER 1 PLUS

BV1-6-10

HEIGHT ADJUSTABLE WORKSTATION WITH COMPUTER MOUNT



Ergonomic



Easy to clean



Small footprint



8 years warranty



Quick setup



Locking wheels

Datasheet:
Click or Scan




FEATURES

**ALL THE FEATURES OF
BEAVER1, PLUS:**

KEYBOARD TRAY AND MOUSE PAD

Comfortable keyboard tray with a left or right hand mouse pad included. Enough space to use the mouse without limitations.



COMPUTER MOUNT WITH POWER STRIP

Safe and secure mount for Tower/Desktop PC. The straps keep the PC from falling.

You can power every device you need with the four-plug power strip.



BEAVER 2

BV2

ALL-IN-ON COMPUTER WORKSTAND



Ergonomic



Easy to clean



Small footprint



8 years warranty



Quick setup



Locking wheels

Datasheet:

Click or Scan



FEATURES

TILT/ROTATE DISPLAY MOUNT

With this solution you can adjust the screen freely. UP to 24" displays.

Tilt: $\pm 15^\circ$
Rotate: $\pm 90^\circ$



SCANNER HOLDER AND HANDLE

We created a solution for moving the stand easily. The scanner holder is on its front, just at the tip of your fingers.



KEYBOARD SHELF WITH PSU HOLDER

We created a big and comfortable shelf to place the keyboard with a mouse. In the PSU holder, you can place power supplies and chargers for your devices or even a 15" laptop.



BEAVER 2 PLUS

BV2-6-4-10

COMPUTER MOUNT WORKSTAND



Ergonomic



Easy to clean



Small footprint



8 years warranty




Quick setup



Locking wheels

Datasheet:
Click or Scan



FEATURES

**ALL THE FEATURES OF
BEAVER2, PLUS:**

COMPUTER MOUNT

Safe and secure mount for Tower/
Desktop PC.
The straps keep the PC from falling.



POWER STRIP

You can power every device
you need with the four-plug
power strip.



BEAVER 3

BV3-6-1

ALL-IN-ONE COMPUTER WORKSTAND



Ergonomic



Easy to clean



Small footprint



8 years warranty




Quick setup



Locking wheels

Datasheet:
Click or Scan



FEATURES

TILT/ROTATE DISPLAY MOUNT


With this solution you can adjust the screen freely. UP to 24" displays.

Tilt: $\pm 15^\circ$
Rotate: $\pm 90^\circ$




SCANNER HOLDER AND HANDLE

We created a solution for moving the stand easily. The scanner holder is on its front, just at the tip of your fingers.



INTERNAL CABLE MANAGEMENT, POWER SUPPLY HOLDER

Hide all cables inside the column body. We added several holes for cables to enter and exit. Storage for the required power supplies.



BEAVER 4

BV4-4-3

LAPTOP CART

Suitable for 13" to 17" laptops.



Ergonomic



Easy to clean



Small footprint



8 years warranty



Quick setup



Locking wheels

Datasheet:

Click or Scan



FEATURES

LEFT-RIGHT HAND MOUSEPAD

Sliding mouse tray can be hidden under the worktop. Comfortable use for both left and right hand.



INTRAORAL SCANNER HOLDER AND SHELF

Small and clever solutions for intraoral scanner. We added a shelf for additional devices.

LAPTOP SECURITY

Optional laptop security to keep the laptop safe from theft and accidental falling.

The protector can be easily mounted on the worktop with knobs.



BEAVER 5

2311-0618

SMALL DEVICE CART



Ergonomic



Easy to clean



Small footprint



8 years warranty



Quick setup



Locking wheels

Datasheet:

Click or Scan



FEATURES

WORKTOP WITH HANDLE

Basic worktop for general use. You can place any instruments on it. We added a handle to make it easy to move.



SHELF

Height adjustable shelf to manage comfortable access. By loosening the screws one can adjust the height of the shelf. Folded up edges prevent objects to fall off.



BEAVER 6

BV6-4-9

PIEZOSURGERY CART



Ergonomic



Easy to clean



Small footprint



8 years warranty




Quick setup



Locking wheels

Datasheet:

Click or Scan



FEATURES

MULTIPLE SHELVES WITH POWER SUPPLY HOLDER

We created large shelves with storage under, where you can place your devices. There's enough space for piezosurgery tools without limitations.



HANDLE

We added a handle to the cart to let it move it easier.

CABLE MANAGEMENT

Multiple holes allow you to run the cables in the column, down to the bottom of the cart.

POWER STRIP

A four-plug powerstrip located on the bottom of the cart to power all you devices needed. Holes added for cable management.





www.dental-carts.com

Medical Cart Company Kft. - 2051 Hungary Biatorbágy, Budai út 16.

Phone: +36 23 920 790 - Email: info@medicalcartcompany.hu

www.medicalcartcompany.hu