

ULTRASOUND CART

SOLUTIONS 2025

Medical Cart
COMPANY



GENERAL FEATURES



EASY TO CLEAN SURFACES

Clean surfaces without dirt traps.
One wipe cleaning.



QUICK SETUP

Available assembled or
flat-packed, easy to assemble.



8 YEARS WARRANTY

Eight-year warranty for
the care-free treatment.



ERGONOMIC

Designed for easy
and comfortable use.



SMALL FOOTPRINT

Developed to fit narrow doors,
ease of movement.



LOCKABLE WHEELS

Simple locking mechanism
operated by foot.



OUR CORE PRODUCTS

CHA4-WS12

BASIC HEIGHT ADJUSTABLE CART

- Height adjustment with gas spring
- Large acryl worktop with handle and holder cuts for Gel bottle and Ultrasonic device
- Display mounting option
- Mounting options on the outer column

Datasheet
Scan the QR



BUG-US0

BASIC STATIC CART

- Tablet or Display mount
- Small acryl worktop with handle and holder cuts for Ultrasonic device
- Mounting options all along the column channel

Datasheet
Scan the QR



3-0050

Lightweight cart with Tablet mount and small acryl worktop.



Datasheet



Scan the QR

3-0072

Height adjustable cart with acryl worktop and Display mounting option.



Datasheet



Scan the QR

3-0240

Height adjustable cart with aluminium worktop and accessories.



Datasheet



Scan the QR

3-0225

Lightweight cart with Tablet arm, acryl worktop, and accessories.



Datasheet



Scan the QR

3-0226


Height adjustable cart with Tablet arm, acryl worktop, and accessories.



Datasheet



Scan the QR



WORKTOP OPTIONS

1-0370

Material:
Acryl

Dimensions:
W: 570mm
D: 475mm
THK: 8mm



1-1111

Material:
Acryl


Dimensions:
W: 350mm
D: 350mm
THK: 8mm



2-0277

Material:
Aluminium


Dimensions:
W: 560mm
D: 485mm
THK: 8mm



WS9

Material:
Acryl

Dimensions:
W: 620mm
D: 560mm
THK: 8mm




WS12

Material:
Acryl

Dimensions:
W: 560mm
D: 485mm
THK: 8mm



**Custom
design
on request**



ACCESSORIES

TABLET ARM



SNA1-WS2-L2B-TM2

Datasheet
Scan the QR



TABLET MOUNT



TM2, TM2-XL

Datasheet
Scan the QR



DISPLAY MOUNT



DM2

Datasheet
Scan the QR



LAPTOP SECURITY



LM2

Datasheet
Scan the QR



GEL HOLDER



GH1, GH3

Datasheet
Scan the QR



SCANNER HOLDER



SR1, SR2, SR3

Datasheet
Scan the QR



PSU HOLDER




PSH 1/2/7/9

Datasheet
Scan the QR



STORAGE



ST2, ST3

Datasheet
Scan the QR



CABLE HOOK



CH1, CH2

Datasheet
Scan the QR



DIN RAIL



DR1

Datasheet
Scan the QR



CUSTOM COLORS

We offer the choice of custom colors for the aluminium components from **Tiger 3D metallic range**.



Visit the Tiger webshop and select the color you desire:
<https://www.tiger-coatings.com>



Medical Cart

C O M P A N Y



Medical Cart Company Kft.

2051 Hungary

Biatorbágy, Budai út 16.

Phone: +36 23 920 790

Email: info@medicalcartcompany.hu

www.medicalcartcompany.hu