

# MEDICAL CART

## SOLUTIONS

### 2025

 Medical Cart  
COMPANY



# INTRODUCTION

We care about the planet and we use materials that fit this thought. Jungle is both the symbol of this approach and the expression of the versatility of our products. The cart ranges are named after animals based on their characteristics.

For the 2024 Catalogue please click or scan the QR code!



## GENERAL FEATURES



### ERGONOMIC

Designed for easy and comfortable use.



### QUICK SETUP

Available assembled or flat-packed, easy to assemble.



### 8 YEARS WARRANTY

Eight-year warranty for the care-free treatment.



### EASY TO CLEAN SURFACES

Clean surfaces without dirt traps. One wipe cleaning.



### SMALL FOOTPRINT

Developed to fit narrow doors, ease of movement.




### LOCKABLE WHEELS

Simple locking mechanism operated by foot.



# BUDGIE RANGE


The most lightweight  
and economic member  
of the product family.





# GIRAFFE RANGE

Vertically packed range  
in case there is no need  
for a work surface.





# RHINO RANGE

Our strongest carts  
for heavy duty  
purposes.





# ZEBRA RANGE

Sit to stand work stations  
with electric lift  
and power options.








# CHAMELEON RANGE

Moving workstation  
that is suitable  
to sit to stand use.



# **DISPLAY RANGE**

Carts for display mounting  
up to 65 inches.







# KINGAROO RANGE


Carts of increased stability  
for multiple devices.





# SNAKE RANGE

Strong and flexible arms  
to provide patient care  
without distance limits.





# DOLPHIN RANGE

Water resistant hardware  
for easy disinfection in  
medical environment.



# TAKE A STAND

Display and tablet stands for those who value comfort and customization. Tilt and height adjustable design providing an optimal viewing angle. Available in a wide range of colors to match any environment and style.



# EASY TO USE EASY TO ADAPT

Comfort and functionality hand in hand.  
We know you have changing needs, and we take them seriously. The carts have modular design that adapts to all your requirements.

## YOU DREAM IT WE MAKE IT!

Realise your ideas and you can count on us to produce them for you. Logos, labels, graphics. We make your brand visible.

Check our visual  
customisation  
options.

**Scan the QR code:**





# EXAMPLE OF A CART

Rhino range cart with possible accessories on the left, fully equipped cart on the right.  
For all carts, accessories and features please see the catalogue in the QR code!

## 1 WORK SURFACE

- 540x620 mm
- Acrylic, available in aluminium

## 3 KEYBOARD TRAY

- Wrist support
- Mouse tray + holder
- Can be hidden under the work surface

## 5 BASE

- B1 Base
- W4 100 mm wheels
- Aluminium with black painting

## 2 DISPLAY MOUNT

- Height adjustable
- Up to 27" displays
- Tilt and rotate functions

## 4 DRAWER

- Push to open
- Shelf on top

## 6 TRANSFORMER TRAY

- Potential equalisation pin on housing
- 260x124x305 mm internal dimensions



## 7 DIN RAIL

- 300x10x25 mm

## 8 RUBBER INSERTS

- For scanners

## 9 COLUMN

- C5 Column
- Fixed height
- Aluminum


## 10 COMPUTER MOUNT

- Variable width

## 11 SHELF

- Custom shape
- Custom painting






Display and tablet  
mount options



Power options



Electric and gas  
spring powered  
lifting options




Potential  
equalisation  
solutions




Holders for PC,  
devices  
and equipments



Power supplies,  
battery pack  
and transformers



Work surface options,  
available in colored  
aluminum



Built-in power  
strips available



Scan the QR!





# Medical Cart COMPANY

**Medical Cart Company Kft.**

2051 Hungary

Biatorbágy, Budai út 16.

Phone: +36 23 920 790

Email: [info@medicalcartcompany.hu](mailto:info@medicalcartcompany.hu)

[www.medicalcartcompany.hu](http://www.medicalcartcompany.hu)